

Title:



Absolute Configuration Determination Report

GENERAL INFORMATION	
Customer	BioTools, Inc.
Sales Order Number	
Sample code (Our ref.)	Penicillin-G
Sample description (Your ref.)	Penicillin-G
VCD-spectrometer	ChirallIR w/ DualPEM
Report prepared by	Jordan Nafie
Report validated and signed by	Rina K Dukor
Date	September 9, 2021
RESULTS	
Absolute Configuration of Penicillin G is (3S,5R,6R) Confidence Level: 99%	
MEASUREMENT PARAMETERS	
Concentration	8mg / 125uL
Solvent	d6-DMSO
Resolution	4 cm ⁻¹
PEM setting	1400 cm ⁻¹
Number of scans/Measurement time	8 hours
Sample cell	BaF ₂
Path length	100 μm
CALCULATION DETAILS	
Gaussian version	Gaussian 09
Total low-energy conformers used for Boltzmann sum	5
Methodology and basis set for DFT calculations	B3LYP / 6311Gdp w/ CPCM (DMSO)
Enantiomer used for calculation	3S,5R,6R
Total calculated conformers	9
Number of low-energy conformations shown in report	3
COMMENTS	
<p>The confidence level is a measure of the degree of congruence between a calculated and measured spectrum. If identical spectra are being compared the confidence level is 100%. The confidence level (CL) is not the likelihood that the assignment is correct. Rather it's a measure of quality or degree of agreement between calculated and measured spectra. With a CL of 99% for this molecule, the visual agreement between measured and calculated spectra is very good – this is a high confidence assignment. Since the sample obtained was the potassium salt, this was measured in d6-DMSO, and calculated as the salt form.</p>	

Title:

Absolute Configuration Determination Report

Structure of Penicillin G:



CompareVOA Results.

Please note: In this plot the frequency scaling factor is not applied.

Title:

Absolute Configuration Determination Report

Table 1. Numerical comparison describing the similarity in the range of 1100 - 1825 cm^{-1} between the calculated IR and VCD spectra for the **(3S,5R,6R)** enantiomer at the B3LYP / 6311Gdp w/ CPCM (DMSO) level and the observed IR and VCD spectra for **Penicillin G**.

Cal. (1100-1825 cm^{-1})	Numerical comparison	Observed Penicillin G
(3S,5R,6R)	scaling factor	0.977
	IR similarity (%)	72.2
	^a Σ (%)	75.5567
	^b Δ (%)	71.809
	Confidence Level (%)	99


^a Σ : single VCD similarity, gives the similarity between the calculated and observed VCD spectra.

^b Δ : enantiomeric similarity index, gives the difference between the values of Σ for both enantiomers of a given diastereoisomer.

Title:

Absolute Configuration Determination Report

Penicillin G in d6-DMSO




IR (lower frame) and VCD (upper frame) spectra of **Penicillin G** in d₆-DMSO; 100 μ m path-length cell with BaF₂ windows; 8 h collection for enantiomer and solvent; instrument optimized at 1400 cm⁻¹. Solvent subtracted IR and VCD spectra are shown. Uppermost trace is the VCD noise spectrum.

Title:

Absolute Configuration Determination Report

Penicillin G Measured vs. Calculated (3S,5R,6R)






IR (lower frame) and VCD (upper frame) spectra **observed** for **Penicillin G** (left axes) compared with Boltzmann-averaged spectra of the **calculated** conformations for the **(3S,5R,6R)** configuration, (right axes).

Title:


Absolute Configuration Determination Report

3 LOWEST ENERGY CONFORMERS - (3S,5R,6R) Configuration:



Title:

Absolute Configuration Determination Report



Plot of ESI (similarity of correct enantiomer minus incorrect enantiomer to calculated) vs SNS (overall similarity of correct enantiomer to calculated) for a library of correct assignments verified independently by X-Ray other method (Black X marks). **Red X** is Penicillin G.